CHAPTER V

CONCLUSION

The results of the present investigation can be summarized as follows:

- Nitrate reductase activity level in normal condition (0.5 M NaCl) is higher than in salt stress condition.
- Nitrate reductase activity level in normal condition containing sodium nitrate
 as nitrogen source is higher than in other nitrogen sources: L-glutamine,
 ammonium chloride. Nitrate can function as an efficient inducer for the
 appearance of nitrate reductase activity level.
- 3. L-glutamine has little effect on growth but affects nitrate reductase activity.
- 4. Ammonium chloride inhibits both growth rate and nitrate reductase activity.
- 5. Nitrate reductase from A. halophytica is localized in cytoplasmic fraction.
- 6. Nitrate reductase from A.halophytica can be purified to homogeneity by four steps including ultracentrifugation, 20 60 % ammonium sulfate precipitation, DEAE-Toyopearl chromatography and Bio-Gel hydroxyapatite chromatography with 15.7 % yield and 406 purification folds.
- The enzyme has a relative molecular weight of monomer about 58 kDa as analyzed by SDS – polyacrylamide gel electrophoresis.
- The purified nitrate reductase from A. halophytica has Km value of nitrate about 465 μM and Vmax value is 32 nmol/min/mg protein in methyl viologen assay.

- The purified nitrate reductase from A. halophytica can used ferredoxin as physiological electron donor.
- 10. The inhibition by various agents reveals that p-Chloromercuribenzoate, iodoacetamide, N-Ethylmaleimide, and KCN can inhibit more than 80 % nitrate reductase activity.
- 11. In the presence of 40 mM ClO₃, the enzyme activity is about 50 % loss.
- 12. Addition 40 mM of nitrite, the product of nitrate reductase activity, has no effect on nitrate reductase activity.