Good and Stenhouse, 22)1966 reported an improved bioassay for TSH by the application of a design balance for residual effects. The modified method has a better increased sensitivity and precision.

MATERIALS

Chemicals, biologicals and animal

- 1. Albino mice* of both sexes, about 10 weeks of age at a weight of 20-25 g.
- 2. <u>Serum samples</u> of normal subjects, myxedema and thyrotoxicosis patients were obtained from the Division of Nuclear Medicine, Department of Radio-logy, Siriraj Hospital.
- 3. Iodine-125 (2 mCi/ml of sodium iodide in M/10 aqueous sodium thiosulphate, from the Radiochemical Centre, Amersham, England) was diluted to 500 uCi/lom
- 4. Thytropar (thyrotropin) obtained from Armour Pharmaceutical Company, Chicago, Illinois 60690, U.S.A., was freshly prepared before use.
- 5. <u>L-Thyroxine</u> (3,3',5,5'-tetraiodo-L-thyronine) obtained from Sigma Chemical Company, St. Louis, U.S.A., stored below 0°C, was diluted to 0.5 ug/ml and freshly prepared before use.
- 6. 3.3'.5-triiodo-L-thyronine (T₃) obtained from Sigma Chemical Company, St. Louis, U.S.A., stored below 0°C, added to the drinking water in the concentration of 0.5 ug/ml and freshly prepared before use.

^{*}Kindly provide by Prof. Dr. Sukhum Patrakom, M.R.V. Dr. Puthipongse Voravuthi, Dr. Pachara Visutakul and Dr. Siri Swasdhikosol.

7. Dog biscuit* of low iodine content (100 ug/Kg) was obtained from W.Menz & Co. Ltd., Galway Avenue, Marsden, South Australia.

The following iodine-free diet may be used in place of the dog biscuit.

Ingredients of Leblond: **

Maize meal	70	%	=	46	lbs.		
Wheat gluten	16'	%	=	102	lbs.		
Yeast	10	%	=	6	lbs.	12	oz.
CaCO ₃	1	%	=	284	g		
NaCl	1	%	=	284	g		

Equipments and glasswares

- 1. Test tubes 10x75 mm and test tube racks
- 2. Heparinized capillary tubes
- 3. Surgical blades No. 24
- 4. Spot test plates
- 5. pipettes 0.1 ml
- 6. Microsyringes (Terumo) 100 ul
- 7. Disposable syringes 1.0 ml
- 8. Electronic Desk-Top Computer (Olivetti Programma 101)
 Cod. Mark.: 02-02-03/9-0-E, Cod. Mag.: 3940758 V
- 9. Automatic Gamma Well Counter System (Nuclear Chicago Model 4230 and 4233).



^{*} Kindly provided by Prof. Basil Hetzel of Monash University, Australia.

^{**}Quoted from Radioisotope Unit, Auckland Hospital, New Zealand.