

REFERENCES

- Jacobsen, S., Fritz, H.G., Degee, Ph., Dubois, Ph., and Jerome, R. (2000) Single-step reactive extrusion of PPLA in a corotating twin-screw extruder promoted by 2-ethylhexanoic acid tin(II) salt and triphenylphosphine. *Polymer*, 41, 3395-3403.
- Jacobsen, S., Fritz, H.G., Degee, Ph., Dubois, Ph., and Jerome, R. (2000) New developments on the ring opening polymerization of polylactide. *Industrial Crops and Products*, 11, 265-275.
- Raquez, J.M., Degee, Ph., Nabar, Y., Narayan, R., and Dubois, Ph. (2006) Biodegradable materials by reactive extrusion: from catalyzed polymerization to functionalization and blend compatibilization. *Comptes Rendus Chimie*, 9, 1370-1379.
- Schmack, G., Jehnichen, D., Vogel, R., Tandler, B., Beyreuther, R., Jacobsen, S., and Fritz, H.G. (2001) Biodegradable fibres spun from poly(lactide) generated by reactive extrusion. *Journal of Biotechnology*, 86, 151-160.
- Hyon, S.H., Jamshidi, K., and Ikada, Y. (1997) Synthesis of polylactides with different molecular weights. *Biomaterials*, 18, 1503-1508.
- Janata, M., Masar, B., Toman, L., Vlcek, P., Latalova, P., Brus, J., and Holler, P. (2003) Synthesis of novel types of graft copolymers by a “grafting-from” method using ring-opening polymerization of lactones and lactides. *Reactive & Functional Polymers*, 57, 137-146.
- Jin, H.J., Chin, I.J., Kim, M.N., Kim, S.H., and Yoon, J.S. (2000) Blending of poly(L-lactic acid) with poly(cis-1,4-isoprene). *European Polymer Journal*, 36, 165-169.
- Odian, G. (2004) *Principles of Polymerization*. New Jersey: John Wiley & Sons.
- Shi, D., Yang, J., Yao, Z., Wang, Y., Huang, H., Jing, W., Yin, J., and Costa, G. (2001) Functionalization of isotactic polypropylene with maleic anhydride by reactive extrusion: mechanism of melt grafting. *Polymer*, 42, 5549-5557.
- Rosales, C., Marquez, L., Perera, R., and Rojas, H. (2003) Comparative analysis of reactive extrusion of LDPE and LLDPE. *European Polymer Journal*, 39, 1899-1915.

- Zhu, L., Tang, G., Shi, Q., Cai, C., and Yin, J. (2006) Neodymium oxide co-catalyzed melt free radical grafting of maleic anhydride onto copolypropylene by reactive extrusion. *Reactive & Functional Polymers*, 66, 984-992.
- Eguiburu, J.L., Fernandez-Berridi, M.J., and Roman, J.S. (1996) Graft copolymers for biomedical applications prepared by free radical polymerization of poly(L-lactide) macromonomers with vinyl and acrylic monomers. *Polymer*, 37(16), 3615-3622.
- Vergnes, B. and Berzin, F. (2006) Modeling of reactive systems in twin-screw extrusion: challenges and applications. *Comptes Rendus Chimie*, 9, 1409-1418.
- Gopakumar, T.G., Ponrathnam, S., Lele, A., Rajan, C.R., and Fradet, A. (1999) In situ compatibilisation of poly(phenylene sulphide)/wholly aromatic thermotropic liquid crystalline polymer blends by reactive extrusion: morphology, thermal and mechanical properties. *Polymer*, 40, 357-364.
- Ryner, M., Stridsberg, K., albertsson, A.C., Schenck, H., and Svensson, M. (2001) Mechanism of ring-opening polymerization of 1,5-dioxepan-2-one and L-lactide with stannous 2-ethylhexanoate: a theory study. *Macromolecule*, 34, 3877-3881.
- Piao, L., deng, M., Chen, X., Jiang, L., and Jing, X. (2003) Ring-opening polymerization of ϵ -caprolactone and L-lactide using organic amino calcium catalyst. *Polymer*, 44, 2331-2336.
- Kelar, K. and Jurkowski, B. (2000) Preparation of functionalized low-density polyethylene by reactive extrusion and its blends with polyamide 6. *Polymer*, 41, 1055-1062.
- Cartier, H. and Hu, G.H. (2001) A novel reactive extrusion process for compatibilizing immiscible polymer blends. *Polymer*, 42, 8807-8816.

CURRICULUM VITAE

Name: Ms. Patcharakamon Nooeaid

Date of Birth: October 15, 1984

Nationality: Thai

University Education:

2002-2005 Bachelor Degree of Chemistry, Faculty of Science, Mahidol University, Bangkok, Thailand

Presentations:

1. Nooeaid, P. and Magaraphan, R. (2007, December 4-7) Catalytic Extrusion of Polylactide/Ethylene Vinyl Alcohol Bioplastic Film. Poster presented at The 10th Polymer Pacific Conference, Kobe, Japan.

